

## **FRCPath Immunology Part 2 practical examination**

### **Station 4 – Data Analysis**

Two example questions are given, including the images that would have been provided on a laptop as Supplementary Material, followed by the answers.

*Note 1: Where data is provided that candidates may want to annotate to help calculate their answers, as in question 2 below, the data may be provided on paper to assist with this.*

*Note 2: Candidates are usually expected to do calculations and so should bring a basic calculator to the exam.*

## Example question 1

Your lab has been alerted by the manufacturer of your rheumatoid factor assay that they will be ceasing manufacture of this assay as a consequence of the new In Vitro Diagnostics Regulations (IVDR). Your lab has an alternative platform that can run this assay. Figure 1 (A,B) shows verification and comparison data for this assay.

Your laboratory's acceptance criteria for assay verification are:

Intra-assay precision <5%

Inter-assay precision <10%

### Figure 1A: Precision data

Table 1: Manufacturer's data for precision

	Mean (IU/ml)	Intra-assay precision	Inter-assay precision
Control 1	27.9	2.7	5.4
Control 2	69.3	2.2	5.2
Control 3	582.7	5.1	4.8


Table 2: Laboratory's verification data for intra-assay precision using the assay QCs

	Low QC	Mid QC	High QC
Mean (n=10)	18.15	26.13	33.57
sd	0.34	0.70	0.49
%CV	1.88	2.68	1.46

Table 3: Laboratory's verification data for inter-assay precision using the assay QCs

	Low QC	Mid QC	High QC
Mean (n=10)	17.33	24.21	31.63
sd	0.90	1.23	1.47
%CV	5.19	5.08	4.72

**Figure 1B: Comparison data**



\*%difference between methods =  
 $(\text{alternative assay result} - \text{current assay result}) / \text{current assay result} \times 100$   
 Blue dotted line is line of best fit.  
 Red dotted line is mean % difference

- (a) Based on the data shown in Figure 1A, would you consider the precision data for this assay acceptable and why (3 marks)
- (b) What do you understand by the term verification? (1 mark)
- (c) Based on the data shown in Figure 1B, comment on the findings of this comparison study (4 marks)
- (d) Reflecting on the information in Figures 1A and 1B, and using your knowledge of this analyte, comment on the QCs used in this assay (3 marks)
- (e) Give one impact of the IVDR on UK labs (1 mark)

## Example question 2

Your laboratory would like to reduce the use of radioactivity in the department. You are comparing the clinical performance of an alternative ELISA-based method (Assay 1) with your existing radioimmunoassay (Assay 2) for TSH receptor antibodies. Review the data in Figure 2 and then answer the questions below:

### Figure 2:

#### A. Assay information

*Assay 1 is an ELISA for TSH receptor antibodies*

*Assay 2 is a radioimmunoassay for TSH receptor antibodies*

#### B. Reference ranges:

	Assay 1	Assay 2
Negative	<0.4 U/L	<1.0 U/L
Borderline	0.4-1.0 U/L	1.0-1.5 U/L
Positive	>1.0 U/L	>1.5U/L

#### C. Manufacturer's information:

*Assay 1: Clinical sensitivity = 95% (n=108), Clinical specificity = 100% (n=137)*

*Assay 2: Clinical sensitivity = 92% (n=50), Clinical specificity = 100% (n=242)*

#### D. Laboratory Comparison data

No.	Assay 1 (U/L)	Assay 2 (U/L)	Assay 1 interpretation	Assay 2 interpretation	Diagnosis of Graves disease (Yes / No)
1	<0.40	<1.0	NEG	NEG	No
2	<0.40	<1.0	NEG	NEG	No
3	3.92	3.8	POS	POS	Yes
4	>30	>40	POS	POS	Yes
5	<0.40	<1.0	NEG	NEG	No
6	6.08	5.3	POS	POS	Yes
7	1.35	2.7	POS	POS	Yes
8	5.08	4.1	POS	POS	Yes
9	2.22	2.5	POS	POS	Yes
10	4.10	3.1	POS	POS	Yes
11	1.82	2.1	POS	POS	No
12	<0.40	<1.0	NEG	NEG	No
13	2.96	2.3	POS	POS	Yes
14	2.23	<1.0	POS	NEG	No
15	4.62	3.5	POS	POS	Yes
16	2.83	2	POS	POS	Yes
17	4.70	3.4	POS	POS	Yes
18	6.39	5.9	POS	POS	Yes
19	2.70	1.8	POS	POS	Yes
20	<0.40	<1.0	NEG	NEG	No
21	22.67	24.2	POS	POS	Yes
22	<0.40	<1.0	NEG	NEG	No
23	<0.40	<1.0	NEG	NEG	No
24	2.01	<1.0	POS	NEG	No
25	5.85	4	POS	POS	Yes
26	<0.40	<1.0	NEG	NEG	No
27	14.58	13.7	POS	POS	Yes
28	<0.40	<1.0	NEG	NEG	No

- Comment on the qualitative agreement between the two assays (2 marks)
- Calculate the clinical sensitivity and specificity for Graves disease for the two assays, stating the equations used and showing your working (6 marks)
- Calculate the positive predictive value of Assay 1 and Assay 2 (2 marks)
- Comment on the potential clinical consequence of the difference in positive predictive value between the two assays (1 mark)

- e) What biochemistry test could you look at to assess or 'sense-check' TSH receptor antibody results, and what result would you expect in a patient who is strongly positive for TSH receptor antibodies for the first time? (2 marks)

## ANSWERS:

### Question 1

- a) Yes
  - Meet manufacturer's stated precision (1 mark)
  - Meets lab's acceptable precision criteria (1 mark)
- b) Verification CONFIRMS a manufacturer's stated performance characteristics (1 mark)
- c) Good linearity (1 mark)
  - Alternative assay has a negative bias compared to current assay/ -11.7% (1 mark)
  - Worse at lower concentrations (1 mark)
  - One outlier, or similar (1 mark)
- d) Low medium and high QCs all at similar concentration (1 mark)
  - QCs should cover dynamic range of assay (1 mark)
  - One QC should be close to clinical cutoff (1 mark)
- e) One of: (1 mark)
  - Harder to do in-house assays
  - Loss of repertoire - companies withdrawing kits
  - Must use commercial assays where available

### Question 2

- a) 26/28 samples / 93% agreement / 2 discrepant samples (1 mark)
  - Reasonable agreement / Would need to investigate clinical case for 2 discrepant samples (1 mark)
- b) Sensitivity equation =  $\text{True pos} / (\text{True pos} + \text{false neg}) \times 100$  (1 mark)
  - Specificity equation =  $\text{True neg} / (\text{True neg} + \text{false pos}) \times 100$  (1 mark)
  - Assay 1 sensitivity =  $16 / (16 + 0) \times 100 = 100\%$  (1 mark)
  - Assay 1 specificity =  $9 / (9 + 3) \times 100 = 75\%$  (1 mark)
  - Assay 2 sensitivity =  $16 / (16 + 0) \times 100 = 100\%$  (1 mark)
  - Assay 2 specificity =  $11 / (11 + 1) \times 100 = 92\%$  (1 mark)
- c) Assay 1 =  $16 / (16 + 3) \times 100 = 84\%$  (1 mark)
  - Assay 2 =  $16 / (16 + 1) \times 100 = 94\%$  (1 mark)

- d) [The lower PPV of assay 1 (ELISA) / more false positives] could lead to unnecessary investigation for Graves disease / incorrect diagnosis of Graves disease (1 mark)
- e) TSH concentration / thyroid function tests (1 mark)  
Expect to see low TSH if TSH receptor antibodies are positive (1 mark)