



UK Health
Security
Agency

UK Standards for Microbiology Investigations

Identification of *Actinomyces* and reclassified species



National Institute for Health and Care Excellence (NICE) has renewed accreditation of the process used by the UK Health Security Agency to produce UK Standards for Microbiology Investigations (UK SMIs). The renewed accreditation is valid until 30 June 2026 and applies to guidance produced using the processes described in 'UK Standards for Microbiology Investigations Development Process' (2021). The original accreditation term began on 1 July 2011.

Consultation between 7 September 2023 to 22 September 2023

Acknowledgments

UK Standards for Microbiology Investigations (UK SMIs) are developed under the auspices of UKHSA working in partnership with the partner organisations whose logos are displayed below and listed on [the UK SMI website](#). UK SMIs are developed, reviewed and revised by various working groups which are overseen by a steering committee (see [the Steering Committee page on GOV.UK](#)).

The contributions of many individuals in clinical, specialist and reference laboratories who have provided information and comments during the development of this document are acknowledged. We are grateful to the medical editors for editing the medical content.

UK SMIs are produced in association with:



Displayed logos correct as of March 2021

UKHSA publications gateway number:

Contents

| | |
|--|----|
| Acknowledgments | 2 |
| Contents | 3 |
| Amendment table | 4 |
| 1 General information | 5 |
| 2 Scientific information..... | 5 |
| 3 Scope of document | 5 |
| 4 Introduction | 5 |
| 5 Technical information and limitations..... | 6 |
| 6 Safety considerations | 6 |
| 7 Target organisms | 7 |
| 8 Identification | 8 |
| 9 Reporting | 14 |
| 10 Referral to reference laboratories..... | 15 |
| Algorithm: Identification of <i>Actinomyces</i> and reclassified species | 16 |
| References..... | 17 |

Consultation between 7 September 2023 to 22 September 2023

Amendment table

Each UK SMI document has an individual record of amendments. The amendments are listed on this page. The amendment history is available from standards@ukhsa.gov.uk.

Any alterations to this document should be controlled in accordance with the local document control process.

| | |
|-------------------------------------|---|
| Amendment number/date | x/dd.mm.yy |
| Issue number discarded | |
| Insert issue number | |
| Anticipated next review date | dd.mm.yy |
| Section(s) involved | Amendment |
| Title | The title has been changed from 'Identification of Anaerobic <i>Actinomyces</i> species' to 'Identification of <i>Actinomyces</i> and reclassified species' |
| | |
| | |

Consultation between 7 September 2023 to 22 September 2023

1 General information

[View general information](#) related to UK SMIs.

2 Scientific information

[View scientific information](#) related to UK SMIs.

3 Scope of document

This UK Standards for Microbiology Investigations (UK SMI) document describes the identification of *Actinomyces* species and includes routine culture, MALDI-TOF MS and gene sequencing for identification. This document does not provide information on antimicrobial susceptibility testing of *Actinomyces* species.

Some of the *Actinomyces* species have been reclassified, and the updated nomenclature of these species have been included in this document for reference.

UK SMIs should be used in conjunction with other relevant UK SMIs.

4 Introduction

4.1 Taxonomy and characteristics

Actinomyces species are gram-positive, filamentous, microaerophilic to facultative anaerobes with high G-C DNA content. The genus *Actinomyces* is in the family *Actinomycetaceae* of the order *Actinomycetales* which belongs to the phylum *Actinomycetota* (*Actinobacteria*) - one of the largest and most diverse phyla among the bacteria (1-3).

Currently, there are 36 *Actinomyces* species validly published with the correct nomenclature and taxonomic status (2). Please refer to the most up to date nomenclature available as changes within the genus are commonplace and may occur following the publication of this UK SMI.

The genus was revised in 2018 and a number of clinically significant species have been designated to alternative genera (see Table 1). For simplicity the original nomenclature has been retained throughout this UK SMI.

Actinomyces species are opportunistic pathogens that form part of the usual microbiota of humans, they typically colonise the oral cavity, gastrointestinal tract, and female urogenital tract. Actinomycosis is a relatively rare and generally polymicrobial infection caused by *Actinomyces* species especially in immunocompromised individuals (4,5).

Identification of *Actinomyces* and reclassified species

The main causative agent is *Actinomyces israelii* but other *Actinomyces* species have also been reported including *Actinomyces odontolyticus*, *Actinomyces viscosus*, *Actinomyces gerencseriae* and *Actinomyces naeslundii*. The clinical presentation of the actinomycosis can vary depending on the severity and site of infection (6-11).

Species within this group have also been identified as clinically significant pathogens in breast abscesses and other non-classical actinomycosis infections (12), of particular note *Arachnia propionica* is associated with canaliculitis (13).

5 Technical information and limitations

Advancements in technology and gene sequencing have significantly contributed to the evolving taxonomy of *Actinomyces*, leading to numerous reclassification and changes in nomenclature of these species (2,3,14). Refer to section 7, table 1 for the reclassified *Actinomyces* species to date.

The complex taxonomy of *Actinomyces* species can lead to uncertainty and inconsistency in their identification, therefore it is important that clinicians stay up to date with the latest taxonomic revisions and resources and incorporate them into their interpretation of laboratory results. It is also important that the databases of identification methods such as MALDI-TOF MS and 16S rRNA gene sequencing reflect any changes in the taxonomy of *Actinomyces* (15).

6 Safety considerations

The section covers specific safety considerations (16-37) related to this UK SMI, and should be read in conjunction with the general [safety considerations on GOV.UK](https://www.gov.uk/safety-considerations).

Actinomyces species and associated Gram-positive species are Hazard group 2 organisms. The processing of diagnostic samples should be carried out at Containment Level 2.

Refer to current guidance on the safe handling of all organisms documented in this UK SMI.

Laboratory procedures that give rise to infectious aerosols must be conducted in a microbiological safety cabinet (25).

The above guidance should be supplemented with local COSHH and risk assessments. Compliance with postal and transport regulations is essential.

7 Target organisms

Table 1. *Actinomyces* species reported to have caused human infection (38-54)

| Previous nomenclature | Current nomenclature |
|--|--|
| <i>Actinomyces israelii</i> | <i>Actinomyces israelii</i> |
| <i>Actinomyces graevenitzii</i> | <i>Actinomyces graevenitzii</i> |
| <i>Actinomyces gerencseriae</i> | <i>Actinomyces gerencseriae</i> |
| <i>Actinomyces naeslundii</i> | <i>Actinomyces naeslundii</i> |
| <i>Actinomyces odontolyticus</i> | <i>Schaalia odontolytica</i> |
| <i>Actinomyces viscosus</i> | <i>Actinomyces viscosus</i> |
| <i>Actinomyces funkei</i> | <i>Schaalia funkei</i> |
| <i>Actinomyces europaeus</i> | <i>Gleimia europaea</i> |
| <i>Actinomyces urogenitalis</i> | <i>Actinomyces urogenitalis</i> |
| <i>Actinomyces meyeri</i> | <i>Schaalia meyeri</i> |
| <i>Actinomyces neuii</i> | <i>Winkia neuii</i> |
| <i>Actinomyces neuii</i> subsp. <i>neuii</i> | <i>Winkia neuii</i> subsp. <i>neuii</i> |
| <i>Actinomyces neuii</i> subsp. <i>anitratus</i> | <i>Winkia neuii</i> subsp. <i>anitrata</i> |
| <i>Actinomyces radingae</i> | <i>Schaalia radingae</i> |
| <i>Actinomyces turicensis</i> | <i>Schaalia turicensis</i> |
| <i>Actinomyces radidentis</i> | <i>Actinomyces radidentis</i> |
| <i>Actinomyces cardiffensis</i> | <i>Schaalia cardiffensis</i> |
| <i>Actinomyces oricola</i> | <i>Actinomyces oricola</i> |
| <i>Actinomyces nasicola</i> | <i>Bowdeniella nasicola</i> |
| <i>Actinomyces massiliensis</i> | <i>Actinomyces massiliensis</i> |
| <i>Actinomyces johnsonii</i> | <i>Actinomyces johnsonii</i> |
| <i>Actinomyces dentalis</i> | <i>Actinomyces dentalis</i> |
| <i>Actinomyces hongkongensis</i> | <i>Pauljensenia hongkongensis</i> |
| <i>Actinomyces hominis</i> | <i>Gleimia hominis</i> |
| <i>Actinomyces oris</i> | <i>Actinomyces oris</i> |
| <i>Actinomyces timonensis</i> | <i>Actinomyces timonensis</i> |
| <i>Actinomyces georgiae</i> | <i>Schaalia georgiae</i> |

Table 2. Other organisms which may be misidentified as *Actinomyces* species

| Previous nomenclature | Current nomenclature |
|--|----------------------------|
| <i>(Pseudo)Propionibacterium propionicum</i> | <i>Arachnia propionica</i> |
| <i>Scardovia wiggisiae</i> | <i>Scardovia wiggisiae</i> |
| <i>Nocardia</i> species | <i>Nocardia</i> species |

8 Identification

Culture-based methods remain the gold standard for identification, however laboratories have increasingly integrated faster identification techniques such as MALDI-TOF MS. There is also a growing shift towards molecular methods for identification. However, these techniques require specialised laboratories, trained staff and expensive reagents which may not be available to all routine laboratories.

8.1 Culture methods

Culture methods provide presumptive identification of *Actinomyces* species based on colony morphology, Gram stain and other phenotypic characteristics followed by identification via MALDI-TOF MS and if necessary, confirmation using 16S rRNA gene sequencing.

8.1.1 Bacterial growth media

Some *Actinomyces* species are fastidious and slow growing (for example, *A israelii*, *A. gerencseriae*) and require enriched medium, with growth enhanced by the addition of carbon dioxide. Anaerobic conditions are favoured but some species can be cultured aerobically or in air plus 5 to 10% CO₂ (55). The optimum growth temperature is 35 to 37°C (54,55). Colonies may appear after 3 to 7 days of incubation, but detection may require 10 to 14 days of incubation (54). Refer to section 8.1.3, table 3 for colony morphology.

Note: The majority of *Actinomyces* species are facultative anaerobes, with the exception of *A israelii*, *A gerencseriae* and *A meyeri* which are strict anaerobes (56).

8.1.1.1 Primary isolation media

Fastidious anaerobic agar or equivalent agar **without neomycin** incubated anaerobically at 35 to 37°C for 5 to 10 days (55).

Note: Many *Actinomyces* species may be inhibited by neomycin.

8.1.1.2 Selective media

Actinomyces selective agar with metronidazole 10mg/L and nalidixic acid 30mg/L (deep fill) incubated anaerobically at 35 to 37°C for 5 to 10 days (55). Growth in air and in air plus 5 to 10% CO₂ is variable. Broth enrichment is rarely beneficial.

Note: Some species may require longer incubation

8.1.2 Colonial appearance

Actinomyces species exhibit various appearances (54,57). Only a few species produce the classic breadcrumb/molar tooth colonies, refer to section 8.2, table 3. The majority are white or grey in colour, with some producing pigmentation following prolonged incubation periods.

8.2 Microscopic appearance

Actinomyces colonies can be examined directly and/or following Gram staining to assess their colony and cellular morphology. Members of the genus demonstrate considerable variation in cellular and colony morphology, which can make their recognition and identification challenging (54,57,58). A combination of different laboratory techniques can aid in their accurate identification.

Note: Clinical specimens such as pus can also be examined directly under the microscope to detect sulphur granules suggesting Actinomycosis (55)

8.2.1 Gram stain

Refer to TP 39 - Staining procedures

Branching, beaded, filamentous or diphtheroid-shaped or coccobacillary Gram positive bacilli.

Note: *Propionibacterium* and *Cutibacterium* species are pleomorphic bacilli that may appear to branch. *Mycardia* species are morphologically indistinguishable from *Actinomyces* species on Gram stain (59).

Table 3. Microscopic and colonial morphology of *Actinomyces* species (2)

The information here provides general characteristics of colony appearance, which can vary among different strains and culture conditions.

| Species | Colonies | Comments |
|------------------------|--|--------------|
| <i>A. israelii</i> | White to cream, breadcrumb, or molar tooth, gritty, pitting | Slow growing |
| <i>A. gerensceriae</i> | Bright white, breadcrumb, or molar tooth, pitting and softer than <i>A. israelii</i> | Slow growing |

Identification of *Actinomyces* and reclassified species

| | | |
|---|---|---|
| <i>A. naeslundii</i> | White, cream, or pinkish, smooth, convex, entire edged | Occasional rough forms occur. Acid production may affect viability in older cultures. |
| <i>A. odontolyticus</i> | Cream to red, smooth, convex, entire edged | Old colonies may be dark brown. Acid production may affect viability in older cultures. |
| <i>A. meyeri</i> | Small, white, smooth, convex, entire edged | Slow growing |
| <i>A. georgiae</i> | White or cream, smooth, convex, entire edged | |
| <i>A. neuui</i> sub sp. <i>neuui</i> and <i>anitratus</i> | White or cream, smooth, convex, entire edged | |
| <i>A. radingae</i> | Grey to white, semi-translucent, smooth, low convex, entire edge | |
| <i>A. turicensis</i> | Grey, semi translucent, smooth, low convex, entire edged | |
| <i>A. europaeus</i> | Whitish, semi translucent, smooth, low convex, entire edged | |
| <i>A. graevenitzi</i> * | White pronounced molar tooth or smooth, convex | Red fluorescence. Rough and smooth forms occur together. Old colonies may become dark brown |
| <i>A. radidentis</i> | Cream to pink, smooth, convex, entire edged | Old colonies may become red |
| <i>A. urogenitalis</i> | Cream to-pink, with darker rings, smooth | Old colonies may become red. Acid production may affect viability in older cultures. |
| <i>A. funkei</i> | Grey, semi translucent, opaque centre (fried egg), low convex, entire edged | |
| <i>A. cardiffensis</i> | Cream to pink, smooth, convex, entire edged | |
| <i>A. nasicola</i> | White or grey, smooth, convex, entire edged | |
| <i>A. oricola</i> | White, breadcrumb, pitting on the agar | |
| <i>A. viscosus</i> | The two types of colonies: large and smooth colonies with V, Y and T configurations or small and rough colonies with short branching filaments. | |
| <i>A. johnsonii</i> | Colonies are similar to <i>A. naeslundii</i> | Acid production may affect viability in older cultures. |
| <i>A. oris</i> | Colonies are similar to <i>A. naeslundii</i> | Acid production may affect viability in older cultures. |
| <i>A. massiliensis</i> | white, pinpoint, circular and shiny with entire edges | |

Issued by the Standards Unit, Specialised Microbiology and Laboratories, UKHSA

Identification | ID 15 | Issue number: dh+ | Issue date: dd.mm.yy| Page: 10 of 22

Identification of *Actinomyces* and reclassified species

| | | |
|----------------------------|---|--|
| <i>A. dentalis</i> | tiny, white and breadcrumb-like and pitting on the agar | |
| <i>A. hongkongensis</i> | non-haemolytic, pinpoint colonies | |
| <i>A. hominis</i> | white–greyish, convex, entire edges | |
| <i>A. timonensis</i> | α-haemolytic, pinpoint, circular, white, dry and embedded in the agar | |
| <i>P. propionicum</i> * | Off white to buff, breadcrumb, gritty, pitting, or smooth, convex, entire edged | Red fluorescence, rough and smooth forms occur together |
| <i>Scardovia wiggisiae</i> | White to cream, breadcrumb, or molar tooth, gritty, pitting. | Very slow growing. May not identify by conventional methods. |

*Colonies of *A. graevenitzii* and *P. propionicum* on blood containing media fluoresce red under long-wave (366 nm) UV illumination.

8.3 Matrix-assisted laser desorption/ionisation - time of flight mass spectrometry (MALDI-TOF MS)

MALDI-TOF MS is currently used as the primary method for the identification of *Actinomyces* species, while 16S rRNA gene sequencing serves as an alternative identification approach or is employed as a confirmatory method.

However, the changing taxonomy and reclassification of *Actinomyces* species poses a challenge for species-level identification (2,3,14). The reference database used by MALDI-TOF MS instruments need to be reviewed regularly and updated with the evolving nomenclature and classification of *Actinomyces* species to ensure accurate identification (15,58,60,63).

Also, closely related *Actinomyces* species exhibit similar protein profiles, making it difficult to distinguish between them, which can lead to misidentification (15,60,61). Additional confirmatory testing with molecular methods such as 16S rRNA gene sequencing may be required.

In addition, the performance of MALDI-TOF MS can be influenced by colony age and morphology. Older colonies or colonies with drier and 'chalkier' morphology can negatively impact the quality of the mass spectral output and subsequently reduce the accuracy of identification of *Actinomyces* species (58,64).

Overall, MALDI-TOF MS is a valuable tool for the identification of *Actinomyces* species, but the factors that limit its effectiveness need to be taken into consideration and tackled by the continual improvement of reference databases and sample preparation processes as well as confirmatory testing when required.

8.4 Further identification

8.4.1 Biochemical tests and commercial identification systems

Biochemical tests including commercial identification kits provide basic biochemical information that can aid in the identification of *Actinomyces* species. However, relying solely on these tests is insufficient for accurate identification of *Actinomyces* species. Therefore, these tests are either not considered reliable for the identification of *Actinomyces* species or employed as part of a multi-step approach that combines alternative identification techniques to achieve more accurate results.

Refer to manufacturer's guidance or Manual of Clinical Microbiology for the biochemical properties of *Actinomyces* species and associated gram-positive species (65).

8.4.1 Indole test

Refer to TP 19 – Indole test

Actinomyces species are spot indole negative

Note: *Propionibacterium/Cutibacterium acnes* is indole positive

8.4.2 Catalase test

Refer to TP 8 – Catalase test

All *Actinomyces* species are catalase negative except *Actinomyces viscosus*, *Actinomyces neuii* subsp *neuii*, *Actinomyces neuii* subsp *anitratus* and *Actinomyces radicidentis* and *Actinomyces horvathi* (57,66).

Laboratories should follow manufacturers' instructions and rapid tests and kits should be validated and be shown to be fit for purpose prior to use.

Results should be interpreted with caution and in conjunction with other test results. To achieve accurate results with biochemical tests, it is advisable to use taxonomic keys and not rely on the identification given by the code. This is because the databases contain out of date information, are incomplete or due variation in reaction strengths occasional weak enzymatic and sugar fermentation reactions (67,68). This is particularly true as molecular techniques enable more species to be identified than was previously possible (68).

8.4.2 Molecular methods

Molecular techniques have made identification of many species more rapid and precise than is possible with phenotypic techniques. The routine implementation of molecular methods can be challenging, as not all clinical laboratories have access to the different molecular methods.

Identification of *Actinomyces* and reclassified species

Therefore, in such cases significant isolates identified by MALDI-TOF MS should be sent to appropriate reference labs for further testing and confirmation of results if required.

8.4.2.1 16S rRNA gene sequencing

The utilisation of 16S rRNA gene sequencing has transformed the identification and taxonomy of *Actinomyces* species (3,14).

Although MALDI-TOF MS is the primary technique for the identification of *Actinomyces* species, 16S rRNA gene sequencing is employed in laboratories as a confirmatory method to resolve any discrepancies and provide more accurate identification of *Actinomyces* species.

However, 16S rRNA gene sequencing also has limitations, these include the challenge of differentiating between some closely related *Actinomyces* species such as *Actinomyces naeslundii*, *Actinomyces viscosus* and *Actinomyces oris*, which can lead to misidentification (58). Within 16S sequencing databases erroneous sequences exist that may appear as the top match. Therefore, it is important to corroborate results using a different identification method such as MALDI-TOF MS and to continue improving and expanding the reference databases to include a comprehensive range of *Actinomyces* species including novel species to reflect current classifications (15,58).

The implementation of 16S rRNA gene sequencing for *Actinomyces* identification can be challenging, as not all clinical laboratories have access to this sequencing technique. In this case any significant isolates should be sent to appropriate reference laboratory.

8.4.2.2 Next generation sequencing (NGS)

With the increased availability of NGS technologies, there may be a shift towards their utilisation for the identification and prediction of antimicrobial susceptibility of *Actinomyces* alongside other target pathogens in future. However currently these technologies are largely restricted to reference units.

Whole genome sequencing (WGS) is routinely used by UKHSA and has greatly improved surveillance capabilities and monitoring trends in antimicrobial resistance. WGS has replaced traditional phenotypic and polymerase chain reaction (PCR) methods for routine surveillance. It also has high discriminatory power and can provide in-depth genetic analysis and identification (69). Therefore, it has the potential to be an alternative to techniques like MALDI-TOF MS for the identification of *Actinomyces* species and associated Gram-positive species.

8.5 Storage and referral

If required, inoculate the pure isolate into anaerobic broth culture or scrape pure growth on to a charcoal transport swab (or anaerobic transport swab) for referral to the appropriate reference laboratory.

9 Reporting

9.1 Infection Specialist

Inform the infection specialist of presumptive or confirmed *Actinomyces* when the request bears relevant information.

9.2 Presumptive identification

If appropriate growth characteristics, colonial appearance and Gram stain of the culture are demonstrated and the isolate is metronidazole non-susceptible.

9.3 Confirmation of identification

If MALDI-TOF MS and/or 16S rRNA gene sequencing identify *Actinomyces* spp. For further confirmation and identification please refer to section 10.

9.4 Health Protection Team (HPT)

N/A

9.5 UK Health Security Agency

N/A

9.6 Infection prevention and control team

N/A

10 Referral to reference laboratories

For information on the tests offered, turnaround times, transport procedure and the other requirements of the reference laboratory [see user manuals and request forms](#)

Contact appropriate reference laboratory for information on the tests available, turnaround times, transport procedure and any other requirements for sample submission:

[England](#)

[Wales](#)

[Scotland](#)


[Northern Ireland](#)

Note: In case of sending away to laboratories for processing, ensure that specimen is placed in appropriate package and transported accordingly.

Consultation between 7 September 2023 to 22 September 2023

Algorithm: Identification of *Actinomyces* and reclassified species

An accessible text description of this flowchart is provided with this document



The flowchart is for guidance only.

Issued by the Standards Unit, Specialised Microbiology and Laboratories, UKHSA

Identification | ID 15 | Issue number: dh+ | Issue date: dd.mm.yy| Page: 16 of 22

References

An explanation of the reference assessment used is available in the [scientific information section on the UK SMI website](#).

1. Schoch CL and others. 'NCBI Taxonomy: a comprehensive update on curation, resources and tools'. Database (Oxford) 2020: volume 2020. 10.1093/database/baaa062
2. Parte AC and others. 'List of Prokaryotic names with Standing in Nomenclature (LPSN) moves to the DSMZ'. International Journal of Systematic and Evolutionary Microbiology 2020: volume 70, issue 11, pages 5607-12. <https://doi.org/10.1099/ijsem.0.004332>
3. Barka EA and others. 'Taxonomy, Physiology, and Natural Products of Actinobacteria'. Microbiology and Molecular Biology Reviews 2016: volume 80, issue 1, pages 1-43. doi:10.1128/membr.00019-15
4. Yun SS and others. 'Lung abscess by *Actinomyces odontolyticus* and *Parvimonas micra* co-infection presenting as acute respiratory failure: a case report'. Medicine 2019: volume 98, issue 35
5. Ding J and others. 'Actinomycotic primary canalitis: Predisposing factors, clinical characteristics, and treatment outcomes'. European Journal of Ophthalmology: volume 0, issue 0. 10.1177/11206721231163512
6. Thukral R and others. 'Actinomyces: a deceptive infection of oral cavity'. Journal of Korean Association of Oral Maxillofacial Surgeons 2017: volume 43, issue 4, pages 282-5. 10.5125/jkaoms.2017.43.4.282
7. Elborno D and others. 'Case Report: Pelvic Actinomyces in an Adolescent with Microperforate Hymen'. Journal of Pediatric and Adolescent Gynecology 2016: volume 29, issue 2, pages 199-201. 10.1016/j.jpag.2016.01.079
8. Walther K and others. 'Actinomyces neuui Isolated From a 20-Month-Old Girl With Cervical Lymphadenitis'. Journal of the Pediatric Infectious Diseases Society 2014: volume 4, issue 3, pages e32-e7. 10.1093/jpids/piu096
9. Markey CM, Vestal LE 'Actinomyces meyeri: A Rare Cause of Postsurgical Pelvic Actinomycosis'. Case Reports in Obstetrics and Gynecology 2018: volume 2018, pages 3842048. 10.1155/2018/3842048
10. Crisafulli E and others. 'A pulmonary infection by *Actinomyces odontolyticus* and *Veillonella atypica* in an immunocompetent patient with dental caries'. Respirology Case Reports 2019: volume 7, issue 9, pages e00493. <https://doi.org/10.1002/rcr2.493>
11. la Cerda-Vargas MF and others. '(99m)Tc-UBI 29-41 bone SPECT/CT scan in craniofacial *Actinomyces israelii*: Misdiagnosis of cranial bone tumor - A case report'. Issued by the Standards Unit, Specialised Microbiology and Laboratories, UKHSA

Surgical Neurology International 2020: volume 11, pages 442.
10.25259/sni_684_2020

12. Bing AU and others. 'Actinomyces Species Isolated from Breast Infections'.
Journal of Clinical Microbiology 2015: volume 53, issue 10, pages 3247-55.
10.1128/jcm.01030-15

13. Seal DV and others. 'Lacrimal canaliculitis due to *Arachnia* (*Actinomyces*)
propionica'. British Journal of Ophthalmology 1981: volume 65, issue 1, pages 10-3
10.1136/bjo.65.1.10

14. Nouioui I and others. 'Genome-Based Taxonomic Classification of the Phylum
Actinobacteria'. Frontiers in Microbiology 2018: volume 9.
10.3389/fmicb.2018.02007

15. Veloo ACM and others. 'Validation of MALDI-TOF MS Biotyper database
optimized for anaerobic bacteria: The ENRIA project'. Anaerobe 2018: volume 54,
pages 224-30. <https://doi.org/10.1016/j.anaerobe.2018.03.007>

16. Advisory Committee on Dangerous Pathogens. '[The Approved List of Biological
Agents](#)'. Health and Safety Executive 2021. pages 1-39
<https://www.hse.gov.uk/pubns/misc208.htm> ++

17. British Standards Institution (BSI). '[BS EN 12469 - Biotechnology - performance
criteria for microbiological safety cabinets](#)'. 2000.
<https://webstore.ansi.org/standards/bsi/bsen124692000> ++

18. British Standards Institution (BSI). '[BS 5726:2005 - Microbiological safety cabinets.
Information to be supplied by the purchaser and to the vendor and to the installer, and
siting and use of cabinets. Recommendations and guidance](#)'. 2005. pages 1-14. ++

19. Centers for Disease Control and Prevention. 'Guidelines for Safe Work Practices
in Human and Animal Medical Diagnostic Laboratories'. Morbidity and Mortality
Weekly Report (MMWR) Surveillance Summaries 2012: volume 61, pages 1-102 +

20. Department for Transport and others. '[Transport of infectious substances UN2814,
UN2900 and UN3373 Guidance note number 17/2012 \(revision 7\)](#)'. 2013 ++

21. Department of Health. '[Health Protection Legislation \(England\) Guidance](#)'. pages
1-124 2010 ++

22. Gizzie N, Adukwu E 'Evaluation of Liquid-Based Swab Transport Systems against
the New Approved CLSI M40-A2 Standard'. Journal of Clinical Microbiology 2016:
volume 54, issue 4, pages 1152-6. ++ 10.1128/JCM.03337-15

23. Health and Safety Executive. '[Managing risks and risk assessment at work](#)'.
(accessed 28/07/2021) ++

24. Health and Safety Executive. '[Safe use of pneumatic air tube transport systems for pathology specimens](#)'. 2009 ++
25. Health and Safety Executive. '[Control of Substances Hazardous to Health Regulations](#)'. The Control of Substances Hazardous to Health Regulations 2002 (as amended). Approved Code of Practice and guidance L5 (sixth edition). Health and Safety Executive (HSE) Books. 2013 ++
26. Health and Safety Executive. '[Risk assessment: A brief guide to controlling risks in the workplace](#)'. Health and Safety Executive (HSE). 2014 ++
27. Health and Safety Executive. '[Advisory Committee on Dangerous Pathogens. Management and operation of microbiological containment laboratories](#)'. Health and Safety Executive (HSE). 2019 ++
28. Health Services Advisory Committee. '[Safe working and the prevention of infection in clinical laboratories and similar facilities](#)'. Health and Safety Executive (HSE) Books. H 2003 ++
29. Home Office. '[Public Health Act \(Northern Ireland\)](#)'. 1967 Chapter 36. 1967 ++
30. Home Office. '[Anti-terrorism, Crime and Security Act](#)'. 2001 ++
31. Official Journal of the European Communities. 'Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on in vitro diagnostic medical devices'. 1998. pages 1-37 ++
32. Public Health England. '[Laboratory reporting to Public Health England: a guide for diagnostic laboratories](#)'. Public Health of England (PHE). 2020. pages 1-31 ++
33. Scottish Government. '[Public Health \(Scotland\) Act](#)'. 2008 ++
34. The Royal College of Pathologists. '[The retention and storage of pathological records and specimens \(5th edition\)](#)'. pages 1-59. 2015 ++
35. The Welsh Assembly Government. '[Health Protection Legislation \(Wales\) Guidance](#)'. 2010 ++
36. Virell KL and others. 'Comparison of the Copan eSwab System with an Agar Swab Transport System for Maintenance of Fastidious Anaerobic Bacterium Viability'. *Journal of Clinical Microbiology* 2016: volume 54, issue 5, pages 1364-7. [10.1128/JCM.03246-15](#)
37. World Health Organization. '[Guidance on regulations for the transport of infectious substances 2019-2020](#)'. World Health Organisation (WHO). 2019 ++

38. Hall V 'Actinomyces-gathering evidence of human colonization and infection'. *Anaerobe* 2008: volume 14, pages 1-7
39. Schaal KP and Lee HJ 'Actinomycete infections in humans-a review'. *Gene* 1992: volume 115, pages 201-11
40. Sabbe LJ and others. 'Clinical spectrum of infections due to the newly described *Actinomyces* Species *A. Turicensis*, *A. radingae* and *A. europaeus*'. *Journal of Clinical Microbiology* 1999: volume 37, pages 8-13
41. Hall V 'Anaerobic actinomycetes and related organisms'. In: Hawkey SGA, editor. *Principles and practice of clinical bacteriology*, second edition. 2 edition: John Wiley & Sons Limited, Chichester UK; 2006. pages. 575-86
42. Hall V and others. '*Actinomyces oricola* sp. nov., from a human dental abscess'. *IntJ Syst Evol Microbiol* 2003: volume 53, issue Pt 5, pages 1515-8
43. Hall V and others. '*Actinomyces dentalis* sp. nov., from a human dental abscess'. *International Journal of Systematic and Evolutionary Microbiology* 2005: volume 55, issue Pt 1, pages 427-31. 10.1099/ijs.0.63376-0
44. Renvoise A and others. '*Actinomyces massiliensis* sp. nov., isolated from a patient blood culture'. *International Journal of Systematic and Evolutionary Microbiology* 2009: volume 59, issue Pt 3, pages 540-4. 10.1099/ijs.0.001503-0
45. Zautner AE and others. 'Subcutaneous fistulae in a patient with femoral hypoplasia due to *Actinomyces europaeus* and *Actinomyces turicensis*'. *Infection* 2009: volume 37, issue 3, pages 289-91. 10.1007/s15010-008-7392-9
46. Woo PC and others. '*Actinomyces hongkongensis* sp. nov. a novel *Actinomyces* species isolated from a patient with pelvic actinomycosis'. *Systematic and Applied Microbiology* 2003: volume 26, issue 4, pages 518-22. 10.1078/072320203770865819
47. Funke G and others. '*Actinomyces europaeus* sp. nov., isolated from human clinical specimens'. *International Journal of Systematic Bacteriology* 1997: volume 47, issue 3, pages 687-92
48. Hall V and others. '*Actinomyces cardiffensis* sp. nov. from human clinical sources'. *Journal of Clinical Microbiology* 2002: volume 40, issue 9, pages 3427-31
49. Funke G and others. '*Actinomyces hominis* sp. nov., isolated from a wound swab'. *International Journal of Systematic and Evolutionary Microbiology* 2010: volume 60, issue Pt 7, pages 1678-81. 10.1099/ijs.0.015818-0
50. Hall V and others. '*Actinomyces nasicola* sp. nov., isolated from a human nose'. *International Journal of Systematic and Evolutionary Microbiology* 2003: volume 53, issue Pt 5, pages 1445-8

51. Renvoise A and others. 'Actinomyces timonensis sp. nov., isolated from a human clinical osteo-articular sample'. International Journal of Systematic and Evolutionary Microbiology 2010: volume 60, issue Pt 7, pages 1516-21. 10.1099/ijls.0.012914-0
52. Könönen E, Wade WG 'Actinomyces and related organisms in human infections'. Clinical Microbiology Reviews 2015: volume 28, issue 2, pages 419-42. 10.1128/cmr.00100-14
53. Gajdács M and others. 'Microbiological and Clinical Aspects of Cervicofacial Actinomyces Infections: An Overview'. Dentistry Journal (Basel) 2019: volume 7, issue 3. 10.3390/dj7030085
54. Gajdács M, Urbán E 'The Pathogenic Role of Actinomyces spp. and Related Organisms in Genitourinary Infections: Discoveries in the New, Modern Diagnostic Era'. Antibiotics (Basel) 2020: volume 9, issue 8. 10.3390/antibiotics9080524
55. Valour F and others. 'Actinomycosis: etiology, clinical features, diagnosis, treatment, and management'. Infection and Drug Resistance 2014: volume 7, pages 183-97. 10.2147/idr.S39601
56. Mbogning Fonkou MD and others. 'Noncontiguous finished genome sequences and descriptions of Actinomyces ihuae, Actinomyces bouchesdurhonensis, Actinomyces urinae, Actinomyces marseillensis, Actinomyces mediterranea and Actinomyces oralis sp. nov. identified by culturomics'. New Microbes and New Infections 2018: volume 25, pages 30-44. <https://doi.org/10.1016/j.nmni.2018.06.004>
57. Sarkonen N and others. 'Phenotypic identification of Actinomyces and related species isolated from human sources'. Journal of Clinical Microbiology 2001: volume 39, issue 11, pages 3955-61. 10.1128/jcm.39.11.3955-3961
58. Fong P and others. 'Identification and diversity of Actinomyces species in a clinical microbiology laboratory in the MALDI-TOF MS era'. Anaerobe 2018: volume 54, pages 151-8. <https://doi.org/10.1016/j.anaerobe.2018.09.007>
59. Langoya CC and others. 'Nocardia and Actinomyces'. Medicine 2021: volume 49, issue 12, pages 756-9. <https://doi.org/10.1016/j.mpmed.2021.09.006>
60. Alcalá L and others. 'Identifying Anaerobic Bacteria Using MALDI-TOF Mass Spectrometry: A Four-Year Experience'. Frontiers in Cellular and Infection Microbiology 2021: volume 11. 10.3389/fcimb.2021.521014
61. Li Y and others. 'Application of MALDI-TOF MS to rapid identification of anaerobic bacteria'. BioMed Central (BMC) Infectious Diseases 2019: volume 19, issue 1, pages 941. 10.1186/s12879-019-4584-0
62. Kim D and others. 'Performance evaluation of a new matrix-assisted laser desorption/ionization time-of-flight mass spectrometry, ASTA MicroIDSys system, in

bacterial identification against clinical isolates of anaerobic bacteria'. *Anaerobe* 2020: volume 61, pages 102131. <https://doi.org/10.1016/j.anaerobe.2019.102131>

63. La Scola B and others. 'Burden of emerging anaerobes in the MALDI-TOF and 16S rRNA gene sequencing era'. *Anaerobe* 2011: volume 17, issue 3, pages 106-12. <https://doi.org/10.1016/j.anaerobe.2011.05.010>

64. Cuénod A and others. 'Factors Associated With MALDI-TOF Mass Spectral Quality of Species Identification in Clinical Routine Diagnostics'. *Frontiers in Cellular and Infection Microbiology* 2021: volume 11. 10.3389/fcimb.2021.646648

65. Karen CC and Michael AP 'Manual of Clinical Microbiology'. 2023

66. Jousimies-Somer H, Summanen P 'Recent taxonomic changes and terminology update of clinically significant anaerobic gram-negative bacteria (excluding spirochetes)'. *Clinical Infectious Diseases* 2002: volume 35, issue Suppl 1, pages S17-S21

67. Kerttula AM and others. 'Enzymatic/biochemical analysis of *Actinomyces* with commercial test kits with an emphasis on newly described species'. *Anaerobe* 2005: volume 11, pages 99-108

68. Santala AM and others. 'Evaluation of four commercial test systems for identification of actinomyces and some closely related species'. *Journal of Clinical Microbiology* 2004: volume 24, pages 418-426

69. Kahler CM 'Neisseria species and their complicated relationships with human health'. *Microbiology Australia* 2021: volume 42, issue 2, pages 79-83. <https://doi.org/10.1071/MA21024>